

# Preparation of Brain Calmodulin

## Day 1

### Materials (all solutions at 4°C)

1. 4 mM EDTA, 40 mM Tris-HCl, pH 8.0. Need 4 liters as the 4x base of other buffers.
2. Buffer A: 1 mM EDTA, 1 mM DTT, 10 mM Tris-HCl, pH 7.5. Prepare by diluting buffer 1 and adding fresh DTT. Prepare 3 mls per gram of tissue.
3. Buffer B: 100 mM NaCl, 1 mM EDTA, 1 mM DTT, 10 mM Tris-HCl, pH 7.5. Need 8 liters.
4. PBS Solution A, 2 liters, ice cold.
5. Ultrapure (enzyme grade) ammonium sulfate.
6. 58% saturation ammonium sulfate. Dissolve 6.94 g ammonium sulfate in 20 ml of 2 NH<sub>2</sub>SO<sub>4</sub>. Save any extra for day 3.
7. 1 M Tris base, a small amount for titration. Save any extra for day 3.
8. Scalpels and dissection scissors.
9. GSA rotor and bottles.
10. Prechilled Waring blender or Polytron with a large generator.
11. DE-52 column, 2.5x20 cm.

### Procedure

1. Obtain fresh brain from a local slaughter house, place in a plastic bag and chill in an ice bucket.
2. Remove membrane, meninges and blood, while keeping the tissue on ice or in a cold room. Collect tissue in a chilled plastic pan.
3. Weigh tissue. Two cow brains yield about 400 g.
4. Transfer to a chilled beaker. Rinse tissue twice with ice cold PBS.
5. Homogenize brain tissue in 2 volumes of buffer A. Pulse the Waring blender several times then blend for ~60 sec in a cold room at medium speed, or use a Polytron for 20 sec.

6. Centrifuge homogenate in a GSA rotor at 12,000 rpm, 4°C for 60 min.
7. Collect pink clear supernatant into a graduate cylinder and read volume. Transfer to a beaker, adjust pH to 7.0 with 1 N NH<sub>4</sub>OH.
8. Add, slowly, dry ammonium sulfate to 58% saturation while monitoring the pH. Need 34.7 grams per 100 ml of supernatant. Soak the electrode in high salt after this step.
9. Stir slowly on ice for 15-20 min, then load into GSA centrifuge bottles and let sit on ice for 40 min. Centrifuge in a GSA rotor at 12,000 rpm, 4°C, for 60 min.
10. Collect supernatant in a beaker, adjust pH to 4.1 with 58% saturation ammonium sulfate in 2 N H<sub>2</sub>SO<sub>4</sub>. Soak the electrode in high salt and clean thoroughly.
11. Let proteins precipitate by sitting on ice for 60 min. Centrifuge in a GSA rotor at 12,000 rpm, 4°C, for 90 min.
12. Resuspend pellets in buffer B and rinse bottles with buffer B. The total volume of buffer used should be ~10% of the starting tissue weight.
13. Adjust pH to 7.5 with 1 M Tris base. Many but not all proteins go back into solution.
14. Dialyze overnight against 2 liters of buffer B. Change buffer twice.
15. Equilibrate the DE-52 column with buffer B.

## **Day 2**

### **Materials**

1. Buffer B, 2 liters.
2. Buffer B, but with 0.5 M NaCl. Need 600 ml.
3. G-100 Sephadex column, 2.5x90 cm.

### **Procedure**

1. Collect solution from the dialysis bag. Centrifuge in a 50.2Ti or 60Ti rotor at 40,000 rpm, 4°C for 30 min.
2. Collect supernatant, load to the DE-52 column. Wash the column thoroughly until OD<sub>280</sub>=0.0 (~500 ml).

3. Elute the column with a 1200 ml gradient of 0.1-0.5 M NaCl in buffer B at a flow rate of ~90 ml/hr. Collect 5-10 ml fractions. Set sensitivity of the UV monitor at OD 2.0. Run the column for 8-10 hr.

4. Equilibrate a G-100 Sephadex column with buffer B, overnight.

## **Day 3**

### **Materials**

1. Buffer B, 2 liters.
2. 1 M Tris base, a small volume for titration, from day 1.
3. 58% saturation ammonium sulfate in H<sub>2</sub>SO<sub>4</sub>, from day 1.
4. Ultrapure (enzyme grade) ammonium sulfate.

### **Procedure**

1. Run SDS-PAGE of column fractions. Calmodulin should appear at 200-250 mM NaCl.
2. Pool peak calmodulin fractions and measure the volume. Bring ammonium sulfate to 58% saturation (34.7 g / 100 ml), while monitoring and adjusting the pH with NH<sub>4</sub>OH.
3. Bring the pH to 4.1 with 58% saturation ammonium sulfate in H<sub>2</sub>SO<sub>4</sub>. Soak the pH electrode in high salt.
4. Let the solution sit on ice for 10-15 min and centrifuge in a SS34 rotor at 17,500 rpm, 4°C for 60 min.
5. Resuspend/dissolve the pellets in a small volume (e.g. 4 ml for 400 g brain) of buffer B. Adjust the pH to 7.5 with 1 M Tris base. Soak the electrode in high salt and clean thoroughly.
6. Centrifuge the solution in a 50Ti rotor at 40,000 rpm, 4°C for 30 min.
7. Load the solution into the G-100 column. Elute at ~40 ml/hr, collect 5 ml fractions. Set sensitivity of the UV monitor at 0.5.

## **Day 4 on**

### **Materials**

1. Buffer C: 150 mM NaCl, 1 mM CaCl<sub>2</sub>, 1 mM DTT, 10 mM Tris-HCl, pH 7.5. Need 2 liters.

2. Buffer C with 0.3 M NaCl, 200 ml.
3. DE-52 column, 1.5x12 cm.
4. 58% saturation ammonium sulfate in H<sub>2</sub>SO<sub>4</sub>, from day 1.
5. Ultrapure (enzyme grade) ammonium sulfate.
6. 10 mM Tris-HCl, pH 7.5, 2 liter.

## Procedure

1. Run SDS-PAGE of column fractions. Pool peak fractions.
2. Bring the solution to 2 mM CaCl<sub>2</sub> and 150 mM NaCl.
3. Dialyze against buffer C for 6 hr to overnight.
4. Equilibrate the DE-52 column with buffer C.
5. Load the solution into the DE-52 column, wash extensively with buffer C.
6. Elute with a 400 ml gradient of 0.15-0.30 M NaCl in buffer C. Collect 3 ml fractions. Set the sensitivity of UV monitor at  $\leq$  OD 0.2.
7. Run SDS-PAGE of column fractions. Pool peak fractions. It may be useful to collect also less pure fractions.
8. Concentrate calmodulin with 58% saturation ammonium sulfate and pH 4.1 as before.
9. Resuspend calmodulin in a small volume of 10 mM Tris-HCl, pH 7.5.
10. Dialyze overnight against 10 mM Tris-HCl. Change buffer twice.
11. Clarify in a 50Ti rotor at 40,000 rpm, 4°C for 30 min.
12. Drop freeze in liquid nitrogen. Stable for months at -80°C.

## Reference

**W.H. Burgess, D.K. Jemiolo and R.H. Kretsinger** (1980) Interaction of calcium and calmodulin in the presence of sodium dodecyl sulfate. *Biochim. Biophys. Acta.* 623:257-270.

---